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

1. Smith G, Greenwood M, Searle R. Ward nurses' use of wound dressings before and after a bespoke educational programme. Journal of Wound Care Vol. 19., no. 9, Sept 2010; 2. Moore Z, Dowsett C, Smith G, et al. TIME CDST: an updated tool to address the current challenges in wound care. Journal of Wound Care, vol 28, no 3, March 2019: 154-161.

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
\*Exsukatmenge für NPWT geeignete Wunden. †SECURA Range includes SECURA Spray, SECURA Protective wipes, SECURA Foam Applicator | 1. Schultz GS, Sibbald RG, Falanga V, et al. Wound bed preparation: a systematic approach to wound management. Wound Rep Reg (2003);11:1-28. 2. Leaper DJ, Schultz G, Carville K, Fletcher J, Swanson T, Drake R. Extending the TIME concept: what have we learned in the past 10 years? Int Wound J 2012; 9 (Suppl. 2):1-19.

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1. Schultz GS, Sibbald RG, Falanga V et al. Wound bed preparation: a systematic approach to wound management. Wound Rep Reg (2003);11:1-28.

2. Leaper DJ, Schultz G, Carville K, Fletcher J, Swanson T, Drake R. Extending the TIME concept: what have we learned in the past 10 years? Int Wound J 2012; 9 (Suppl. 2):1-19.

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1. Tallis A, et al. Clinical and Economic Assessment of Diabetic Foot Ulcer Debridement with Collagenase: Results of a Randomized Controlled Study. *Clinical Therapeutics*/Volume 35, 11 (2013). 2. Milne CT, et al. A Comparison of Collagenase to Hydrogel Dressings in Wound Debridement. *Wounds* 22 (11):270–274 (2010). 3. Atkin L. A guide to wound debridement. <http://www.independentnurse.co.uk/clinical-article/a-guide-to-wound-debridement/114500> (2016).

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

1. Schultz GS, Sibbald RG, Falanga V et al. Wound bed preparation: a systematic approach to wound management. Wound Rep Reg (2003);11:1-28. 2. Leaper DJ, Schultz G, Carville K, Fletcher J, Swanson T, Drake R. Extending the TIME concept: what have we learned in the past 10 years? Int Wound J 2012; 9 (Suppl. 2):1-19.

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1. Guest, JF. Costs and Consequences of Managing Infection Associated with Venous Leg Ulcers in the UK. Presentation. CATALYST Health Economics Consultants Ltd Northwood, Middlesex, UK. in (2017). 2. Guest JF, Fuller GW, Vowden P. Diabetic foot ulcer management in clinical practice in the UK: costs and outcomes. *Int Wound J* 15, 43–52 (2018). 3. Malone M, et al. The prevalence of biofilms in chronic wounds: a systematic review and meta-analysis of published data. *J Wound Care* 26, 20–25 (2017). 4. Schultz G, et al. Consensus guidelines for the identification and treatment of biofilms in chronic nonhealing wounds. *Wound Repair Regen.* 25, 744–757 (2017). 5. Schierle CF, De la Garza M, Mustoe TA, Galiano RD. Staphylococcal biofilms impair wound healing by delaying reepithelialization in a murine cutaneous wound model. *Wound Repair Regen.* 17, 354–9 (2009). 6. Roche ED, et al. Increasing the presence of biofilm and healing delay in a porcine model of MRSA-infected wounds. *Wound Repair Regen.* 20, 537–43 (2012). 7. Zhao G, et al. Time course study of delayed wound healing in a biofilm-challenged diabetic mouse model. *Wound Repair Regen* 20, 342–352 (2012). 8. Costerton JW, Stewart PS, Greenberg EP. Bacterial biofilms: a common cause of persistent infections. *Science* 284, 1318–22 (1999). 9. Stewart PS, Costerton JW. Antibiotic resistance of bacteria in biofilms. *Lancet* 358, 135–8 (2001). 10. Stewart PS, Rayner J, Roe F, Rees WM. Biofilm penetration and disinfection efficacy of alkaline hypochlorite and chlorosulfamates. *J. Appl. Microbiol.* 91, 525–32 (2001). 11. Phillips PL, et al. Antimicrobial dressing efficacy against mature *Pseudomonas aeruginosa* biofilm on porcine skin explants. *Int. Wound J.* 12, 469–83 (2013). 12. Bjarnsholt T, et al. Why chronic wounds will not heal: a novel hypothesis. *Wound Repair Regen.* 16, 2–10 (2008). 13. Upton D, Solowiej K. Pain and stress as contributors to delayed wound healing. *Wound Pract. Res.* 18, 114–122 (2010).

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
1. Schultz G, et al. Consensus guidelines for the identification and treatment of biofilms in chronic nonhealing wounds. *Wound Repair Regen.* 25, 744–757 (2017). 2. Guest JF, Fuller GW, Vowden P. Venous leg ulcer management in clinical practice in the UK: costs and outcomes. *Int Wound J.* 15, 29–37 (2018). 3. Dowsett C, Bielby A, Searle R. Reconciling increasing wound care demands with available resources. *J. Wound Care* 23, 552, 554, 556–8 passim (2014). 4. Newton H. Reducing MRSA bacteraemias associated with wounds. *Wounds UK* 6, 56–65 (2010). 5. Searle R, Bielby A. Dressing strategies for the management of infected wounds in community wound care : impacts and implications. *Wounds UK* (2010).

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
Idealfeuchtes Wundmilieu

1. Schultz GS, Sibbald RG, Falanga V, et al. Wound bed preparation: a systematic approach to wound management. Wound Rep Reg (2003);11:1-28. 2. Leaper DJ, Schultz G, Carville K, Fletcher J, Swanson T, Drake R. Extending the TIME concept: what have we learned in the past 10 years? Int Wound J 2012; 9 (Suppl. 2):1-19.

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---

1. Tickle J, 2015. Wound exudate assessment and management: a challenge for clinicians. British Journal Of Nursing 24, pp. S38-S43. 2. O'Keefe M. Evaluation of a community based wound care programme in an urban area. Poster presented at EWMA Conference 2006. 3. College of Nursing; Nursing Shortage Preventing Access to Crucial Training. Says RCN. Healthwatch Trafford. 2015. <https://healthwatchtrafford.co.uk/news/nursingshortage-preventing-access-to-crucial-training-says-rcn/>. 4. Maranozov R, Huxley C, Manzoni C, Pike G. 2017. Royal college of nursing employment survey, Institute of employment studies.

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
Clinical learnings adapted from: Dowsett C, et al. Use of PICO to improve clinical and economic outcomes in hard-to-heal wounds. Wounds International (2017) Vol 8 (No. 2): 53-58.

The products used in the T.I.M.E. clinical decision support tool may vary in different markets. Not all products referred to may be approved for use or available in all markets. Please consult your local Smith+Nephew representative for further details on products available in your market. Intended for healthcare professionals outside of the US only.

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
1. Schultz G, Sibbald G, Falanga V, Ayello E, Dowsett C, Harding K, Romanelli M, Stacey M, Teot L, Vanscheidt W. Wound bed preparation: a systematic approach to wound management. Wound repair and regeneration. Vol. 11, No. 2, Supplement March–April 2003. S1-S28.

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1. Lindholm C, Searle R. Wound management for the 21st century: combining effectiveness and efficiency. *Int Wound J.* 2016;13 (Suppl 2):5–15. 2. Guest JF, Vowden K, Vowden P. The health economic burden that acute and chronic wounds impose on an average clinical commissioning group/health board in the UK. *J Wound Care.* 2017;26:292–303. 3. Ousey K, Stephenson J, Barrett S, et al. Wound care in five English NHS Trusts: results of a survey. *Wounds UK.* 2013;9:20–28. 4. Dowsett C, Bielby A, Searle R. Reconciling increasing wound care demands with available resources. *J Wound Care.* 2014;23:552, 554, 556–8. 5. Leeds Community Healthcare NHS Trust RY6, Community Health Services for adult's quality report. May 2016.

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
1. Lindholm C, Searle RJ. Wound management for the 21st century: combining effectiveness and efficiency. Int Wound J. 2016;13 Suppl 2:5–15.15. 2. Jørgensen SF, Nygaard R, Posnett J. Meeting the challenges of wound care in Danish home care. J Wound Care. 2013;22:540–2, 544–545.

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



\*Exsuaratmenge für NPWT geeignete Wunden. †SECURA Range includes SECURA Spray, SECURA Protective wipes, SECURA Foam Applicator | 1. Schultz GS, Sibbald RG, Falanga V, et al. Wound bed preparation: a systematic approach to wound management. Wound Rep Reg (2003);11:1-28. 2. Leaper DJ, Schultz G, Carville K, Fletcher J, Swanson T, Drake R. Extending the TIME concept: what have we learned in the past 10 years? Int Wound J 2012; 9 (Suppl. 2):1-19.

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

1. Guest JF, et al. Health economic burdens that wounds impose on the National Health Service in the UK. BMJ Open. 2015;5:e009283.

2. Guest JF. Costs and consequences of managing infection associated with venous leg ulcers in the UK. Presented at: CATALYST Health Economics Consultants Ltd Northwood; 2017; Middlesex, UK.

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



1. Guest JF, et al. Health economic burden that wounds impose on the National Health Service in the UK. *BMJ Open*. 2015;5:e009283.  
 2. Malone M, et al. The prevalence of biofilms in chronic wounds: a systematic review and meta-analysis of published data. *J. Wound Care*. 2017;26:20–25.  
 3. Guest JF, Fuller GW, Vowden P. Venous leg ulcer management in clinical practice in the UK: costs and outcomes. *Int. Wound J*. 2018;15:29–37.  
 4. Guest JF, Fuller GW, Vowden P. Diabetic foot ulcer management in clinical practice in the UK: costs and outcomes. *Int. Wound J*. 2018;15:43–52.

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[www.smith-nephew.com](http://www.smith-nephew.com)  
◇ Trademark of Smith+Nephew  
All Trademarks acknowledged  
©Dezember 2019 Smith+Nephew  
AWM-Miscellaneous-22309

**Smith+Nephew**

---

Smith+Nephew does not provide medical advice. The information presented is not, and is not intended to serve as, medical advice. For detailed device information, including indications for use, contraindications, precautions and warnings, please consult the product's Instructions for Use (IFU) prior to use. It is the responsibility of healthcare professionals to determine and utilise the appropriate products and techniques according to their own clinical judgment for each of their patients.

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
2 of 4


1. DaCosta RS, Kulbatski I, Lindvere-Teene L, Starr D, Blackmore K, Silver JI, et al. Point-of-care autofluorescence imaging for real-time sampling and treatment guidance of bioburden in chronic wounds: first-in-human results. PLOS one. 2015 Mar 19;10(3). 2. Rennie MY, Lindvere-Teene L, Tapang K, Linden R, Point-of-care fluorescence imaging predicts the presence of pathogenic bacteria in wounds: a clinical study. Journal of Wound Care 2017;26(8):452-460.

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3 of 4

1. DaCosta RS, Kulbatski I, Lindvere-Teene L, Starr D, Blackmore K, Silver JI, et al. Point-of-care autofluorescence imaging for real-time sampling and treatment guidance of bioburden in chronic wounds: first-in-human results. PLOS one. 2015 Mar 19;10(3).

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1. Hill R, Douglas JJ. Real-time bacterial fluorescence imaging guides antimicrobial stewardship in patients with diverse wounds. Proceedings of the Annual Symposium on Advanced Wound Care (SAWC); 2017 Apr 5-9; San Diego, CA. 2. Landis S. Use of fluorescence imaging in visualizing bacteria in chronic ulcers and traumatic soft tissue damage. Proceedings of the Annual Meeting of the Society of Federal Health Professionals (AMSUS); 2016 Nov 29-Dec 2; National Harbor, MD. 3. Ottolino-Perry K, Chamma E, Blackmore KM, Lindvere-Teene L, Starr D, Tapang K, et al. Improved detection of wound bacteria using fluorescence image guided wound sampling in diabetic foot ulcers. Int Wound J. 2017 Feb 28. doi: 10.1111/iwj.12717.

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1. Colin D, Kurring P, Quinlan D, Yvon C. The clinical investigation of an amorphous hydrogel compared with a dextranomer paste dressing in the management of sloughy pressure ulcers. Proceedings of the 5th European conference on advances in wound management, Harrogate, UK, November 1995. 2. Flanagan M. The efficacy of a hydrogel in the treatment of wounds with Non-viable tissue. Journal of Wound Care, 1995 Vol 4(6), 264 - 267. 3. Spooner R. Managing a patients multiple pressure sores. Journal of Wound Care (1993), 2(3), 139 - 141. 4. Ricci E, et al. INTRASITE Gel: indications and results in the therapy of Decubitus Ulcers. Personal Communication 1993. 5. Thomas S. The role of moist wound healing in the management of meningococcal skin lesions: a case study. World Wide Wounds, June 1999. 6. Young T, Williams C, Benbow M, et al. A study of two hydrogels used in the management of pressure sores. 103-106. 7. Bale S, Banks V, Haglestein S, et al. A comparison of two amorphous hydrogels in the debridement of pressure sores. Journal of Wound Care 1998; 7 (2): 65-68. 8. Coates TJ, et al. The effect of gel burns dressings on skin temperature. Emergency Medicine Journal. 2002; 19:224-225. 9. Thomas S, Rowe H, Keats J, Morgan R. The management of extravasation injury in neonates. World Wide Wounds October 1997. (Also in Thomas S, Rowe H, Keats J, Morgan R. A new approach to the management of extravasation injury in neonates. The Pharmaceutical Journal 1987;239: 584). 10. Ricci E, Cassino R, Secreto P. Treatment of necrotic tissue with INTRASITE gel. 11. Regan M. The use of INTRASITE Gel in healing open sternal wounds. Ostomy/Wound management, Journal for Extended Patient Care Management, Volume 38, No. 3, April 1992. 12. Flanagan M. The use of a Modified Formulation of an Amorphous Hydrogel Containing Cross-linked Carboxy Methyl Cellulose (CMC) Polymer in the Management of Wounds that Present with Non-viable Tissue – A Multi-centre Clinical Trial. 13. Westerhof W, Mekkes JR. A pilot study comparing INTRASITE Gel with saline soaked gauze for debridement, Ref CTR90/08 by D Myers. 14. Colin D, Kurring PA, Quinlan D, Yvon C. Managing sloughy pressure sores. Journal of Wound Care, 1996 Vol 5(10), 444 – 446. 15. Thomas S, Fear M. The efficacy of INTRASITE Gel as a debrider of nonviable tissue from wounds treated in the community. Poster presentation at the European Wound Management Conference, Harrogate 1993. 16. Pastre L, Wetzberger S, Andriessen A. INTRASITE Gel as a desloughing agent in the management of patients with wounds requiring debridement in the community setting. 17. De Vincentiis G, Caracciolo G, Anselmi A. INTRASITE Gel in the management of deep second and third degree burns in children [Poster] 1996. [Poster] 1995. 18. Thomas S, Rowe HN, et al. A new approach to the management of extravasation injury in neonates. The Pharmaceutical Journal (1987) Nov: 584 - 585. 19. Ratliff CR. INTRASITE Gel to promote wound healing from necrosis. The Cancer Bulletin 191; 43 (1): 85-6. 20. Gates JL, Holloway GA. A comparison of wound environments. A moist wound environment dressing system versus the traditional normal saline wet-to-dry dressing. Ostomy/Wound Management 1992; 38 (8): 35-37. 21. Boardman M, Mellor K, Neville B. Treating a patient with a heavily exuding malodorous fungating ulcer. Journal of Wound Care. 1993 1993;2(2): 74-6. 22. Van Look R, Verbist A. The treatment of complicated and infected wounds [Poster] 1993. 36. Riggs R, Bale S. Management of necrotic wounds as a complication of histiocytosis X. Journal of Wound Care. 1993 1993;2(5):260-1. 23. Mitchell E. Keeping leg ulcers moist. Nurs Times. 1990 9/19/1990;86(38):67-8. 24. Benbow M, Pearce C. The care of an infant with homozygous protein C deficiency. Journal of Wound Care. 1994 1994;3(1):21-4. 25. Hollinworth H. Managing a patient with an infected foot ulcer. Journal of Wound Care. 1993 1993;2(1):22-6. 26. Meier L. Treatment of a patient with an amputation above the left knee [Poster] 1994. 27. Todd M. The use of an amorphous hydrogel in the management of granulating wounds. Poster at Woundcare Society Meeting, Harrogate, UK, November 1994. 28. Ricci E, Barberid S, Ferraris S, Anselmetti G. The treatment of fistulae with INTRASITE Gel; Interim results [Poster] 1993. 29. Cassino R, Ricci E, Carusone A. Management of heavily exuding wounds requiring debridement [Poster] 1997. 30. Gethin G, Cowman S. Manuka Honey vs. Hydrogel—a prospective, open label, multicenter, randomized controlled trial to compare desloughing efficacy and healing outcomes in venous ulcers. J Clin Nurs. 2009 2/2009;18(3):466-74.

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
1. Granick MS et al, "Efficacy and cost-effectiveness of a high-powered parallel waterjet for wound debridement", Wound Repair & Regeneration, 2006, 14, 394-397. 2. Granick MS et al, "Toward a common language: surgical wound bed preparation and debridement", Wound Repair & Regeneration, 2006, 14, S1-S10. 3. Cubison CS et al, "Dermal preservation using the Versajet hydrosurgery system for debridement of pediatric burns", Burns, 2006, 32, 714-720. 4. Chammas MF, Gurunluoglu R, Carlsen SN, et al. Surgical debridement of mineral pitch and nonviable penile tissue using water-jet power: a preliminary report. BJU International 2008;103(7):974-976. 5. Vanwijck R, Kaba L, Boland S, et al. Immediate skin grafting of sub-acute and chronic wounds debrided by hydrosurgery. Journal of Plastic, Reconstructive & Aesthetic Surgery. 2010;63:544-549.

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
1 of 4

1. Skog E, Arnesjo B, Troeng T, et al. A randomised trial comparing cadexomer iodine and standard treatment in the outpatient management of chronic venous ulcers. *Br J Dermatology* (1983) 109, p77-83. 2. Sundberg J, Meller R. A retrospective review of the use of cadexomer iodine in the treatment of chronic wounds. *Wounds* (1997) 9(3), p68-86. 3. Troeng T, Skog E, Arnesjo B, et al. In Cadexomer Iodine, Fox JA, Fischer H (eds) (1983) p43-50. 4. Holloway GA, Johansen KH, Barnes RW, et al. Multicenter trial of cadexomer iodine to treat venous stasis ulcer. *West J Med* (1989, Jul) 151, p35-38. 5. Ormiston MC, Seymour MTJ, Venn GE, et al. A randomised comparison of cadexomer iodine and a standard treatment in out-patients with chronic leg ulcers. In Cadexomer Iodine, Fox JA, Fischer H (eds) (1983) p63-69. 6. Hansson C. The effects of cadexomer iodine paste in the treatment of venous leg ulcers compared with hydrocolloid dressing and paraffin gauze dressing. *Int J Dermatology* (1998) 37, p390-396. 7. Moberg S, Hoffman L, Grenner M, et al. A randomised trial of cadexomer iodine in decubitus ulcers. *J Am Geriatrics Soc* (1983) 31, p462-465. 8. Lindsay G, Latta D, Lyons KGB, et al. A study in general practice of the efficacy of cadexomer iodine in venous leg ulcers treated on alternate days. *Acta Therapeutica* (1986) 12, p141-148. 9. Johnson A. A combative healer with no ill effect. IODOSORB® in the treatment of infected wounds. *Professional Nurse* (Oct 1991) p60-64. 18. Smith & Nephew Wound Management, Data on file report 0810016, Antimicrobial activity of ACTICOAT FLEX against a broad spectrum of wound pathogens, Driffeld K, October 2008. 19. Smith & Nephew Wound Management, Data on file report 0810012, Antimicrobial activity of ACTICOAT FLEX 7 against a broad spectrum of wound pathogens, Driffeld K, October 2008.

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


1. Wright JB, et al. (1998) Wound management in an era of increasing bacterial antibiotic resistance: A role for topical silver treatment, American Journal of Infection Control: 26(6):572-577. 2. Driffield K. Antimicrobial activity of ACTICOAT Flex 7 against a broad spectrum of wound pathogens. Data on file 0810012, Smith+Nephew. 3. Upton D. (2010) Pain and stress as contributors to delayed wound healing. Wound Practice and Research. 5. Gago M, et al. (2008) A comparison of Three Silver – containing Dressings in the Treatment of Infected, Chronic Wounds 20(10): 273-278. 6. Daubney L. Silver Release Testing of ACTICOAT Flex 7 Dressings. Report reference DS/08/062/R2.

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

1. Wright JB, et al. (1998) Wound management in an era of increasing bacterial antibiotic resistance: A role for topical silver treatment, American Journal of Infection Control: 26(6):572-577. 2. Wright JB, et al. (1998) Efficacy of topical silver against fungal burn wound pathogens, American Journal of Infection Control (1999); 27(4): 344-350. 3. Driffield K. ACTICOAT Flex 3 has antimicrobial activity in 30 minutes. Data on file 0810018, Smith+Nephew. 4. Driffield K. ACTICOAT Flex 7 has antimicrobial activity in 30 minutes. Data on file 0810014, Smith+Nephew. 5. Gago M, et al. (2008) A comparison of Three Silver – containing Dressings in the Treatment of Infected, Chronic Wounds 20(10): 273-278. 6. Westaim Report Ref. #010322. Seven Day Efficacy of ACTICOAT 7 Dressings Against Multiple Organisms. 7. Driffield K. Antimicrobial activity of ACTICOAT Flex 7 against a broad spectrum of wound pathogens. Data on file 0810012, Smith+Nephew. 8. Driffield K. Antimicrobial activity of ACTICOAT Flex 7 dressings in a seven-day repeat challenge. Data on file 0810013, Smith+Nephew. 9. Westaim Report Ref: # 971030. The Antimicrobial Activity of Westaim's ACTICOAT Silver Coated Dressing against Clinically Relevant Organisms over an Extended Period of Time. 10. Driffield K. Antimicrobial Activity of ACTICOAT Flex 3 against a Broad Spectrum of Wound Pathogens. Data on File reference 0810016. 11. Driffield K. Antimicrobial activity of ACTICOAT Flex 3 in a three-day repeat challenge. Data on file 0810017, Smith+Nephew. 12. Upton D. (2010) Pain and stress as contributors to delayed wound healing. Wound Practice and Research. 13. Westaim (Sherritt) Report Ref: 93/001. Broad Spectrum Efficacy. 14. Daubney, L. Silver Release Testing of ACTICOAT Flex 7 Dressings. Report reference DS/08/062/R2. 15. Daubney, L. Silver Release Testing of ACTICOAT Flex 3 Dressings. Report reference DS/08/078/R2.

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
1 of 4

1. Fitzgerald DJ, Renick PJ, Forrest EC, et al. Cadexomer iodine provides superior efficacy against bacterial wound biofilms in vitro and in vivo. *Wound Repair Regen.* 2017;25(1):13-24. 2. Phillips PL, Yang Q, Davis S, et al. Antimicrobial dressing efficacy against mature *Pseudomonas aeruginosa* biofilm on porcine skin explants. *Int Wound J.* 2015;12(4):469-483. 3. Schultz G, Qingpin Y. Microbicidal Effects of Three Daily Treatments of a Carboxymethylcellulose Silver Dressing or a Cadexomer Iodine on Mature Bactericidal Biofilms Grown on Pig Skin Explants. 2016; Florence. 4. Roche ED, Woodmansey EJ, Yang Q, et al. Cadexomer iodine effectively reduces bacterial biofilm in porcine wounds ex vivo and in vivo. *Int Wound J.* 2019;1-10 [in press]. 5. Akiyama H, Oono T, Saito M, Iwatsuki K. Assessment of Cadexomer iodine against *Staphylococcus aureus* Biofilm In vivo and In vitro Using Confocal Laser Scanning Microscopy. *The Journal of Dermatology.* 2004;31:529-534. 6. Fitzgerald DJ, Renick PJ, Forrest EC, et al. Cadexomer iodine provides superior efficacy against bacterial wound biofilms in vitro and in vivo. *Wound Repair Regen.* 2017;25(1):13-24. 7. Phillips PL, Yang Q, Davis S, et al. Antimicrobial dressing efficacy against mature *Pseudomonas aeruginosa* biofilm on porcine skin explants. *Int Wound J.* 2015;12(4):469-483. 8. Oates JL, Phillips CD, Wolcott R, Woodmansey E. Effect of a Cadexomer Iodine Wound Dressing on a Chronic Wound Multi-Species Biofilm Model with Comparison to a Silver Hydrofiber Antibiofilm Dressing. 2016. 9. Hill KE, Malic S, Mckee R, et al. An in vitro model of chronic wound biofilms to test wound dressings and assess antimicrobial susceptibilities. *J Antimicrob Chemother.* 2010;65(6):1195-1206. 10. Kiamco MM, Atci E, Mohamed A, Call DR, Beyenal H. Hyperosmotic Agents and Antibiotics Affect Dissolved Oxygen and pH Concentration Gradients in *Staphylococcus aureus* Biofilms. *Appl Environ Microbiol.* 2017;83(6). 11. O'Meara S, Al-Kurdi D, Ologun Y, Ovington LG, Martyn-St James M, Richardson R. Antibiotics and antiseptics for venous leg ulcers (Review). *Cochrane Database of Systematic Reviews.* 2014(Issue 1 Art. No.: CD003557). 12. O'Meara S, Al-Kurdi D, Ologun Y, Ovington LG. Antibiotics and antiseptics for venous leg ulcers (Review). *Cochrane Database of Systematic Reviews.* 2010(Issue 1 Art. No.: CD003557).

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



1. Roche ED, Woodmansey EJ, Yang Q, et al. Cadexomer iodine effectively reduces bacterial biofilm in porcine wounds ex vivo and in vivo. *Int Wound J.* 2019;1-10 [in press]. 2. Akiyama H, Oono T, Saito M, Iwatsuki K. Assessment of Cadexomer Iodine against Staphylococcus aureus Biofilm In vivo and In vitro Using Confocal Laser Scanning Microscopy. *The Journal of Dermatology.* 2004;31:529-534. 3. Fitzgerald DJ, Renick PJ, Forrest EC, et al. Cadexomer iodine provides superior efficacy against bacterial wound biofilms in vitro and in vivo. *Wound Repair Regen.* 2017;25(1):13-24. 4. Phillips PL, Yang Q, Davis S, et al. Antimicrobial dressing efficacy against mature Pseudomonas aeruginosa biofilm on porcine skin explants. *Int Wound J.* 2015;12(4):469-483. 5. Oates JL, Phillips CD, Wolcott R, Woodmansey E. Effect of a Cadexomer Iodine Wound Dressing on a Chronic Wound Multi-Species Biofilm Model with Comparison to a Silver Hydrofiber Antibiofilm Dressing. 2016. 6. Hill KE, Malic S, Mckee R, et al. An in vitro model of chronic wound biofilms to test wound dressings and assess antimicrobial susceptibilities. *J Antimicrob Chemother.* 2010;65(6):1195-1206. 7. Kiamco MM, Atci E, Mohamed A, Call DR, Beyenal H. Hyperosmotic Agents and Antibiotics Affect Dissolved Oxygen and pH Concentration Gradients in Staphylococcus aureus Biofilms. *Appl Environ Microbiol.* 2017;83(6). 8. Ormiston MC, Seymour MTJ, Venn GE, Cohen RI, Fox JA. Controlled trial of Iodosorb<sup>®</sup> in chronic venous ulcers. *British Medical Journal.* 1985;291:308-310. 9. Hansson C, Persson L-M, Stenquist B, et al. The effects of cadexomer iodine paste in the treatment of venous leg ulcers compared with hydrocolloid dressing and paraffin gauze dressing. *International Journal of Dermatology.* 1998;37(390-396). 10. Troeng T, Skog E, Arnesjo B, et al. A Randomised Multicentre Trial to Compare the Efficacy of Cadexomer Iodine and Standard Treatment in the Management of Chronic Venous Ulcers in Out-Patients. Stuttgart: Schattauer Verlag; 1983. 11. Skog E, Arnesjo B, Troeng T, et al. A randomized trial comparing cadexomer iodine and standard treatment in the out-patient management of chronic venous ulcers. *British Journal of Dermatology.* 1983;109:77-83. 12. Smith+Nephew data on file 2009; #091101. 13. Gottardi W. Disinfection, sterilization and preservation (ed. Block, S. S.) 152-166 (Lea and Febiger, 1991). 14. Cooper RA. Iodine revisited. *Int Wound J.* 4, 124-37 (2007). 15. Johnson A. A combative healer with no ill-effect Iodosorb in the treatment of infected wounds. *Prof. Nurse* 7, 60, 62, 64 (1991).

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<sup>‡</sup> Staphylococcus aureus mature biofilms; <sup>§</sup> MRSA biofilms; <sup>§</sup> Mixed bacterial cultures. Pseudomonas aeruginosa PA01, Staphylococcus aureus Mu50, and Enterococcus faecalis V583; <sup>¥</sup> Treatment every 24 h for 48 h total. 1. Phillips, et al. Int W J. 12: 469-83 (2013) 2. Fitzgerald DJ, et al. Wound Repair Regen. 25: 13-24 (2016); 3. Oates JL, et al. Poster presented at SAWC, Atlanta.(2016); 4. Schultz G, Yang. Poster presented at WUWHS Florence (2016) N.B. Values in red indicate the number of times more effective IODOSORB was demonstrated to be compared to Aquacel Ag+ in that test model.

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1. Smith+Nephew Wound Management laboratory report, reference DS/16/341/R. Analytical and Physical Properties of DURAFIBER Ag Recipe E. Giaffar L, October 2016. 2. Smith+Nephew Wound Management laboratory report, reference DS/15/269/R. DURAFIBER Ag (Recipe E) Testing. Dowler A, December 2015. 3. Forlee M, Rossington A, Searle R. A prospective, open, multicentre study to evaluate a new gelling fibre dressing containing silver in the management of venous leg ulcers. Int. Wound J. 2014;11(4):438-45. 4. Smith+Nephew Wound Management laboratory report, reference DS/12/158/R2. Sequestration Testing of DURAFIBER Ag. Lidster C, April 2014. 5. Smith+Nephew Wound Management data on file report, reference 1412017. Visual demonstration of sequestration of Pseudomonas aeruginosa and Staphylococcus aureus in a DURAFIBER Ag dressing using confocal microscopy. Lumb H, December 2014. 6. Smith+Nephew Wound Management data on file report, reference 1510009. Antimicrobial activity of DURAFIBER Ag against bacteria, yeast and fungi commonly found in wounds over a 7 day period. Lumb H, October 2015. 7. Smith+Nephew Wound Management data on file report, reference 1510008. Antimicrobial activity of DURAFIBER Ag against a broad spectrum of wound pathogens. Lumb H, October 2015. 8. Smith+Nephew Wound Management data on file report, reference 15010010. Antimicrobial activity of DURAFIBER Ag against bacteria, yeast and fungi commonly found in wounds over a 4 hour period. Lumb H, October 2015.

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1 of 2

\*Based on proven performance on wound model data. ‡Based on *in vitro* data. §On average. ||In ALNB pilot evaluation with 29 patients vs clinicians' usual dressings.

1. Smith+Nephew data on file DS.16.469.R. ALLEVYN LIFE Non-Bordered Vertical Wound Model. 2. Smith+Nephew R&D Report Version 2 RD/17/097v2 - Technical description of dressing mode of action of ALLEVYN LIFE Non-Bordered. 3. Smith+Nephew data on file DS17.530.R Comparison of fluid absorption under compression between ALLEVYN LIFE Non-Bordered and Mepilex XT 2017. 4. Smith+Nephew data on file DS17.529.R. Comparison of fluid absorption under compression between ALLEVYN LIFE Non-Bordered and Mepilex 2017. 5. Smith+Nephew data on file 171103;10,20-22 Pilot study data summary.

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
\* When used with appropriate secondary retention

1. Smith+Nephew R&D Report Version 2 WMP/13/507 – Technical description of dressing mode of action of ALLEVYN LIFE Non-Bordered. 2. Smith+Nephew data on file DS17.530.R. Comparison of fluid absorption under compression between ALLEVYN LIFE Non-Bordered and Mepilex XT Smith+Nephew data on file DS.16.334 v2 ALLEVYN LIFE Non-Bordered Fluid Uptake Testing. 3. Rossington A, et al. Clinical performance and positive impact on patient wellbeing of ALLEVYN Life. *Wounds UK* 2013;9(4):91-5. 4. Data on file 1611018 (2016). Bacterial barrier testing of ALLEVYN Life Non-Bordered using modified wet/dry conditions. 5. Smith+Nephew data on file DS.16.328.R ALLEVYN Life Non-Bordered Physical Testing. 2017

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
\*The bacterial properties of the dressing will be compromised if it is cut.

1. Rossington A, et al. Clinical performance and positive impact on patient wellbeing of ALLEVYN Life. Wounds UK 2013;9(4):91-5. 2. Joy H, et al. A collaborative project to enhance efficiency through dressing change practice. Journal of Wound Care 2015;24(7):312,314-7. 3. Data on File Report DS/12/127 - Odour reducing properties of ALLEVYN Life. 14 June 2012. 4. Smith & Nephew 20 June 2016. A Randomised Cross-Over Clinical Evaluation to Compare Performance of ALLEVYN Life and Mepilex® Border Dressings on Patient Wellbeing-Related Endpoints. Internal Report. CE/047/ALF

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\*Wie in-vitro nachgewiesen. References: 1. Smith+Nephew 2019. New ALLEVYN Gentle with Silicone Wound Contact Layer – Physical Properties. Internal Report. RD/19/083/V2 2. Smith+Nephew 2019. Project Peridot (ALLEVYN Gentle V2 with Alston Silicone WCL) Non-Specification Physical Testing for Verification. Internal Report.DS/19/058/R. 3. Smith+Nephew 2018.Volunteer Evaluation of AWM Product Features. Internal Report. EO.AWM.PCSgen.005.v1. 4. Smith+Nephew 2019.Bacterial Barrier testing of Peridot ALLEVYN Gentle dressing verification batches under wet/wet conditions. Internal Report. MB-TM001-29-R. 5. Smith+Nephew 2019.ALLEVYN Gentle 2 vs ALLEVYN Gentle – horizontalwound model using viscous SWF at a flow rate modelling that of a moderately exuding wound. Internal Report. DS/19/054/R. 6. Rafter L, Reynolds T, Rafter M. An audit of patient outcomes in the management of skin tears using silicone dressings. Wounds UK. 2016;12(2):70-78. 7. Carter J. The use of a silicone adhesive dressing (ALLEVYN Gentle Border) on a patient with extremely fragile skin. Paper presented at: Wounds UK; 2008;Harrogate, UK. 8. Smith+Nephew 2014. An Open, Prospective, Randomised Volunteer Trial to Assess the Retention Qualities of 4 Non-bordered Silicone Foam Dressings under Secondary Retention. Internal Report. VTR/HVT079/ST965. 9. Vowden K. An open, prospective randomised, multi-centre clinical evaluation of a hydrocellular silicone foam dressing\* in the management of exuding chronic and acute wounds. Paper presented at: EWMA; 2011; Brussels. 10. Hurd T, Gregory L, Jones A, Brown S. A multi-centre in-market evaluation of ALLEVYN Gentle Border. Wounds UK. 2009;5(3):32-44. 11. Benbow M. Exploring the concept of moist wound healing and its application in practice. British Journal of Nursing. 2008;17 (Tissue Viability Supplement)(15):S4 -S16. 12. Smith+Nephew 2019.Project Peridot (ALLEVYN Gentle V2 with Alston Silicone WCL) Non-Specification Physical Testing for Verification. Internal Report.DS/19/058/R. 13. Smith+Nephew 2019. ALLEVYN Gentle 2 – 7 day horizontal wound model testing using horse serum at a flow rate modelling that of a highly exuding wound. Internal Report. DS/19/051/R. 14. Smith+Nephew 2019.ALLEVYN Gentle 2 – 7 day horizontal wound model testing using horse serum at a flow rate modelling that of a moderately exudingwound. Internal Report. DS/19/052/R. 15. Smith+Nephew 2019.Use of Moisture Vapour Permeability\* (MVP) and Moisture Vapour Transmission Rate\* (MVTR) data to support ALLEVYN Gentle (Project Peridot) providing a moist wound environment which is optimal for healing. Internal Report. EO.AWM.PCS263.001.v1.

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
1. Colin D, Kurring P, Quinlan D, Yvon C. The clinical investigation of an amorphous hydrogel compared with a dextranomer paste dressing in the management of sloughy pressure ulcers. Proceedings of the 5th European conference on advances in wound management, Harrogate, UK, November 1995. 2. Flanagan M. The efficacy of a hydrogel in the treatment of wounds with Non-viable tissue. Journal of Wound Care, 1995 Vol 4(6), 264 - 267. 3. Spooner R. Managing a patient's multiple pressure sores. Journal of Wound Care (1993), 2(3), 139 - 141. 4. Ricci E, et al. INTRASITE Gel: indications and results in the therapy of Decubitus Ulcers. Personal Communication 1993. 5. Thomas S. The role of moist wound healing in the management of meningococcal skin lesions: a case study. World Wide Wounds, June 1999. 6. Young T, Williams C, Benbow M, et al. A study of two hydrogels used in the management of pressure sores. 103-106. 7. Bale S, Banks V, Haglestein S, et al. A comparison of two amorphous hydrogels in the debridement of pressure sores. Journal of Wound Care 1998; 7 (2): 65-68. 8. Coates TJ, et al. The effect of gel burns dressings on skin temperature. Emergency Medicine Journal. 2002; 19:224-225. 9. Thomas S, Rowe H, Keats J, Morgan R. The management of extravasation injury in neonates. World Wide Wounds October 1997. (Also in Thomas S, Rowe H, Keats J, Morgan R. A new approach to the management of extravasation injury in neonates. The Pharmaceutical Journal 1987;239: 584). 10. Ricci E, Cassino R, Secreto P. Treatment of necrotic tissue with INTRASITE gel. 11. Regan M. The use of INTRASITE Gel in healing open sternal wounds. Ostomy/Wound management, Journal for Extended Patient Care Management, Volume 38, No. 3, April 1992. 12. Flanagan M. The use of a Modified Formulation of an Amorphous Hydrogel Containing Cross-linked Carboxy Methyl Cellulose (CMC) Polymer in the Management of Wounds that Present with Non-viable Tissue – A Multi-centre Clinical Trial. 13. Westerhof W, Mekkes JR. A pilot study comparing INTRASITE Gel with saline soaked gauze for debridement, Ref CTR90/08 by D Myers. 14. Colin D, Kurring PA, Quinlan D, Yvon C. Managing sloughy pressure sores. Journal of Wound Care, 1996 Vol 5(10), 444 - 446. 15. Thomas S, Fear M. The efficacy of INTRASITE Gel as a debrider of nonviable tissue from wounds treated in the community. Poster presentation at the European Wound Management Conference, Harrogate 1993. 16. Pastre L, Wetzberger S, Andriessen A. INTRASITE Gel as a desloughing agent in the management of patients with wounds requiring debridement in the community setting. 17. De Vincentiis G, Caracciolo G, Anselmi A. INTRASITE Gel in the management of deep second and third degree burns in children [Poster] 1996. [Poster] 1995. 18. Thomas S, Rowe HN, et al. A new approach to the management of extravasation injury in neonates. The Pharmaceutical Journal (1987) Nov: 584 - 585. 19. Ratliff CR. INTRASITE Gel to promote wound healing from necrosis. The Cancer Bulletin 191; 43 (1): 85-6. 20. Gates JL, Holloway GA. A comparison of wound environments. A moist wound environment dressing system versus the traditional normal saline wet-to-dry dressing. Ostomy/Wound Management 1992; 38 (8): 35-37. 21. Boardman M, Mellor K, Neville B. Treating a patient with a heavily exuding malodorous fungating ulcer. Journal of Wound Care. 1993 1993;2(2): 74-6. 22. Van Look R, Verbist A. The treatment of complicated and infected wounds [Poster] 1993. 36. Riggs R, Bale S. Management of necrotic wounds as a complication of histiocytosis X. Journal of Wound Care. 1993 1993;2(5):260-1. 23. Mitchell E. Keeping leg ulcers moist. Nurs Times. 1990 9/19/1990;86(38):67-8. 24. Benbow M, Pearce C. The care of an infant with homozygous protein C deficiency. Journal of Wound Care. 1994 1994;3(1):21-4. 25. Hollinworth H. Managing a patient with an infected foot ulcer. Journal of Wound Care. 1993 1993;2(1):22-6. 26. Meier L. Treatment of a patient with an amputation above the left knee [Poster] 1994. 27. Todd M. The use of an amorphous hydrogel in the management of granulating wounds. Poster at Woundcare Society Meeting, Harrogate, UK, November 1994. 28. Ricci E, Barberid S, Ferraris S, Anselmetti G. The treatment of fistulae with INTRASITE Gel; Interim results [Poster] 1993. 29. Cassino R, Ricci E, Carusone A. Management of heavily exuding wounds requiring debridement [Poster] 1997. 30. Gethin G, Cowman S. Manuka Honey vs. Hydrogel—a prospective, open label, multicenter, randomized controlled trial to compare desloughing efficacy and healing outcomes in venous ulcers. J Clin Nurs. 2009 2/2009;18(3):46-74.

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1. Dowler A. DURAFIBER Dressing physical properties. Laboratory report DS/10/056/R1 05 2010. 2. Dowler A. Testing of absorption of DURAFIBER and Aquacel<sup>®</sup> while the dressing is under compression. Analytical report DS/10/037/R1 03 2010. 3. Dowler A. Report of testing on lateral wicking of DURAFIBER versus Aquacel and KALTOSTAT<sup>®</sup> dressings. Analytical report DS/10/042/R1 04 2010. 4. Demonstrable. 5. Barret S, et al. Case series evaluation: The use of DURAFIBER on exuding wounds. Wounds UK 2012, Vol 8, No 3. 6. Lumb H. Visual demonstration of sequestration of Pseudomonas aeruginosa and Staphylococcus aureus in a DURAFIBER Dressing using confocal microscopy. Data on file report - 1004001 04 2010. 7. Dowler A. Testing performed as part of DURAFIBER photo shoot. October 2011. Laboratory report DS/11/187/DOF 8. Woodmansey E. The potential for proliferation of bacteria in wound exudate. RR-WMPO6290-40-01 04 2010.

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1. Dowsett C, Hampton J, Myers D, Styche T. Use of PICO to improve clinical and economic outcomes in hard-to-heal wounds. Wounds International. 2017;8, p53-58. 2. Hampton J. Providing cost-effective treatment of hard-to-heal wounds in the community through use of NPWT. Br J Community Nurs, 2015. Suppl Community Wound Care:S1, S16-20. (13). 3. Kirsner R, Dove C, Reyzelman A, Vayser D, Jaimes H. Randomized controlled trial on the efficacy and acceptance of a single-use negative pressure wound therapy system versus traditional negative pressure wound therapy in the treatment of lower limb chronic ulcers (VLU and DFU). Poster 8 presented at Wild on Wounds National Wound Conference. September 12-15, 2018, Las Vegas, USA.

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1. Smith+Nephew 2016. Smith+Nephew's RENASYS<sup>◇</sup> Touch Device – User Validation. Internal Report. DS/16/174/R. 2. Forlee M, Richardson J. An interim analysis of device functionality and usability of RENASYS TOUCH - a new portable Negative Pressure Wound Therapy (NPWT) system. Paper presented at Wounds UK 2016, Harrogate, UK. 3. Forlee M, Van Zyl L, Louw H, Nel J, Fourie N, Hartley R. A randomised controlled trial to compare the clinical efficacy and acceptability of adjustable intermittent and continuous Negative Pressure Wound Therapy (NPWT) in a new portable NPWT system. Paper presented at EWMA, 2018, Krakow, Poland. 4. Smith+Nephew 2016. Smith+Nephew's RENASYS Touch Device Features – Software features. Internal Report. DS.16.179.R. 5. Smith+Nephew 2018. Project Colne Verification of RENASYS TOUCH Connect device battery while in use. 21. DS.17.623.R V2.

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\*n = 71; Chronic wounds = 24, Sub-acute = 37 and Acute = 10; 4-week RCT study period

1. Forlee M, Van Zyl L, Louw H, Nel J, Fourie N, Hartley R. A randomised controlled trial to compare the clinical efficacy and acceptability of adjustable intermittent and continuous Negative Pressure Wound Therapy (NPWT) in a new portable NPWT system. Paper presented at EWMA, 2018. Krakow, Poland. 2. Forlee M, Richardson J. An interim analysis of device functionality and usability of RENASYS<sup>◇</sup> TOUCH - a new portable Negative Pressure Wound Therapy (NPWT) system. Paper presented at Wounds UK 2016. Harrogate, UK.

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94%  
RENASYS

90%  
V.A.C

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1. Hurd T, et al. A Retrospective Comparison of the Performance of Two Negative Pressure Wound Therapy Systems in the Management of Wounds of Mixed Etiology. *Advances in Wound Care*, 2017. 6(1): p. 33-37.

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\*Based on 5 out of 9 wounds responding; wound mean duration prior to study 44 week, study size n=9. †A prospective cohort study of 52 wounds, PICO healed wounds n=14

1. Dowsett C, Hampton J, Myers D, Styche T. Use of PICO to improve clinical and economic outcomes in hard-to-heal wounds. Wounds International. 2017;8, p53–58. 2. Hampton J. Providing cost-effective treatment of hard-to-heal wounds in the community through use of NPWT. Br J Community Nurs, 2015. Suppl Community Wound Care:51, S16–20. (13). 3. Kirsner R, Dove C, Reyzelman A, Vayser D, James H. Randomized controlled trial on the efficacy and acceptance of a single-use negative pressure wound therapy system versus traditional negative pressure wound therapy in the treatment of lower limb chronic ulcers (VLU and DFU). Poster 8 presented at Wild on Wounds National Wound Conference. September 12–15, 2018, Las Vegas, USA.

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4 of 4

1. Dowsett C, Huddleston E, Milne J. Reducing the burden of chronic wounds in the community with single-use NPWT Using disposable negative (2015). 2. Roberts S. In-vitro wound model testing of PICO at a moderate exudate flow rate, March 2011. Data on file reference DS/11/037/R2. 3. Roberts S. In-vitro wound model testing of PICO at a low exudate flow rate, April 2011. Data on file reference DS/11/057/R2. 4. Malmjö M, et al. Biological effects of a disposable, canisterless Negative Pressure Wound Therapy system. ePlasty 2014;14:e15. 5. Lumb H. Bacterial barrier testing (wet-wet) of PICO Dressing with a 7 day test duration against S. marcescens. February 2011 (In-vitro). Data on file report 1102010. 6. Young S, Hampton S, Martin R. Non invasive assessment of Negative Pressure Wound Therapy using high frequency diagnostic ultrasound: oedema reduction and new tissue accumulation (2013) IWJ 10, p383-388. 7. Canonico et al. Therapeutic possibilities of portable NPWT. Initial multidisciplinary observation with the negative pressure therapy device. ACTA VULNOL 2012;10:57-66. 8. Edwards J. Use of a portable Negative Pressure Wound Therapy System (NPWT) (PICO) for split thickness skin grafts; Nov 2012.

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1 of 4

n = 71; Chronic wounds = 24, Sub-acute = 37 and Acute = 10; 4-week RCT study period

1. Smith+Nephew 2016. Smith+Nephew's RENASYS Touch Device – User Validation. Internal Report. DS/16/174/R. 2. Forlee M, Richardson J. An interim analysis of device functionality and usability of RENASYS TOUCH - a new portable Negative Pressure Wound Therapy (NPWT) system. Paper presented at Wounds UK 2016, Harrogate, UK. 3. Forlee M, Van Zyl L, Louw H, Nel J, Fourie N, Hartley R. A randomised controlled trial to compare the clinical efficacy and acceptability of adjustable intermittent and continuous Negative Pressure Wound Therapy (NPWT) in a new portable NPWT system. Paper presented at EWMA, 2018, Krakow, Poland. 4. Smith+Nephew 2016. Smith+Nephew's RENASYS Touch Device Features – Software features. Internal Report. DS.16.179.R. 5. Smith+Nephew 2018. Project Colne Verification of RENASYS TOUCH Connect device battery while in use. 21. DS.17.623.R V2.

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